

ATCC® Human Site-Specific Microbiome Reference Standards



Anna McCluskey, BS, Briana Benton, BS, Juan Lopera, PhD, and Dev Mittar, PhD | ATCC, Manassas, VA 20110

Background

A predominant limitation in microbiome research is the lack of appropriate and relevant standards to control the technical biases introduced throughout the metagenomics workflow. To address this, ATCC has developed a set of genomic DNA and whole cell mock microbial communities from fully sequenced and characterized ATCC strains that represent species found in the oral, skin, gut, or vaginal microbiome. To further enhance the use of these standards and eliminate the bias associated with data analysis, we have also collaborated with One Codex to develop data analysis modules that provide simple output in the form of true-positive, relative abundance, and false-negative scores for 16S rRNA community profiling and shotgun metagenomics sequencing.

Relevance of Strains in Site-specific Microbiome Standards

Table 1. The strains selected for the site-specific microbiome standards were chosen on the basis of their relevance in the normal and typical flora of the oral, skin, gut, and vaginal microbiomes

Species	Gram Stain	Genome Size (Mb)	% GC	Relevance
Oral Microbiome Standard (Genomic DNA: ATCC® MSA-1004™ and Whole Cell: ATCC® MSA-2004™)				
<i>Actinomyces odontolyticus</i>	Positive	2.39	65.5	Anaerobic bacterium associated with dental plaques
<i>Fusobacterium nucleatum</i>	Negative	2.17	27.2	Periodontal pathogen
<i>Haemophilus parainfluenzae</i>	Negative	2.12	39.3	Common oral commensal
<i>Prevotella melaninogenica</i>	Negative	3.17	35.1	Associated with dental caries
<i>Streptococcus mitis</i>	Positive	1.83	40.5	Associated with dental caries
<i>Veillonella parvula</i>	Negative	2.16	38.6	Prevalent on all oral surfaces
Skin Microbiome Standard (Genomic DNA: ATCC® MSA-1005™ and Whole Cell: ATCC® MSA-2005™)				
<i>Acinetobacter johnsonii</i>	Negative	3.47	41.9	Frequently encountered in the skin microbiota
<i>Corynebacterium striatum</i>	Positive	2.72	59.7	Nosocomial opportunistic pathogen
<i>Cutibacterium acnes</i>	Positive	2.49	60.1	Causative agent of acne
<i>Micrococcus luteus</i>	Positive	2.50	73.0	High GC content
<i>Staphylococcus epidermidis</i>	Positive	2.50	32.0	Common cause of nosocomial infections
<i>Streptococcus mitis</i>	Positive	1.83	40.5	Commensal but occasionally pathogenic
Gut Microbiome Standard (Genomic DNA: ATCC® MSA-1006™ and Whole Cell: ATCC® MSA-2006™)				
<i>Bacteroides fragilis</i>	Negative	5.21	43.3	Symbiotic but occasionally opportunistic in the peritoneal cavity
<i>Bacteroides vulgatus</i>	Negative	5.16	42.2	Most common fecal isolate from humans
<i>Bifidobacterium adolescentis</i>	Positive	2.09	59.2	Found in breast-fed newborns
<i>Clostridioides difficile</i>	Positive	4.11	28.6	May colonize the gut following antibiotic therapy
<i>Enterobacter cloacae</i>	Negative	5.31	55.1	Opportunistic pathogen following antibiotic exposure
<i>Enterococcus faecalis</i>	Positive	3.36	37.4	Opportunistic pathogen, produces cytolyisin toxin
<i>Escherichia coli</i>	Negative	4.64	50.6	Typically found in the lower intestine of humans
<i>Fusobacterium nucleatum</i>	Negative	2.17	27.0	Belongs to normal microflora of oral and gastrointestinal tracts
<i>Helicobacter pylori</i>	Negative	1.67	38.9	Associated with peptic ulcers and chronic gastritis
<i>Lactobacillus plantarum</i>	Positive	3.31	44.5	Commonly found in probiotics to regulate intestinal microflora
<i>Salmonella enterica</i>	Negative	4.59	52.1	Not considered to be normal microflora
<i>Yersinia enterocolitica</i>	Negative	4.55	47.0	Foodborne and waterborne pathogen that causes gastroenteritis
Vaginal Microbiome Standard (Genomic DNA: ATCC® MSA-1007™ and Whole Cell: ATCC® MSA-2007™)				
<i>Gardnerella vaginalis</i>	Negative	1.67	41.9	Sexually transmitted pathogen associated with bacteremia and UTIs
<i>Lactobacillus gasseri</i>	Positive	1.89	34.9	Found in the oral, intestinal, and vaginal microflora of humans
<i>Lactobacillus jensenii</i>	Positive	1.67	34.4	Found in the normal flora of the human female urogenital tract
<i>Mycoplasma hominis</i>	Negative	0.67	27.0	Causative agent of urogenital infections
<i>Prevotella bivia</i>	Negative	2.52	39.9	Associated with endometritis and septic arthritis
<i>Streptococcus agalactiae</i>	Positive	2.16	35.4	Associated with septicemia, meningitis, and pneumonia in newborns

Analysis of Site-specific Microbiome Communities




Figure 1. Genomic DNA mock microbial communities can be analyzed via 16S rRNA and shotgun metagenomics sequencing methods. Genomic DNA mixes representing the (A) oral (ATCC® MSA-1004™), (B) skin (ATCC® MSA-1005™), (C) gut (ATCC® MSA-1006™), and (D) vaginal (ATCC® MSA-1007™) microbiomes were analyzed via shotgun metagenomics and 16S rRNA sequencing. Data analyses were performed on the One Codex platform.

DNA Extraction Can Affect Microbiome Analysis Data




Figure 2. Whole cell mock microbial communities can be used to evaluate the impact of DNA extraction on the microbiome analysis data. The Oral Microbiome Genomic DNA Mix (ATCC® MSA-1004™) and genomic DNA extracted from the Oral Microbiome Whole Cell Mix (ATCC® MSA-2004™) were analyzed via whole-genome sequencing and the resulting data were compared. The Oral Microbiome Genomic DNA Mix produced consistent results with the expected abundance of 16.7%. In contrast, the Oral Microbiome Whole Cell Mix produced variable results; this variability may be attributed to differences in DNA extraction efficiency associated with cell wall properties.

Shotgun Metagenomics versus 16S rRNA Analysis




Figure 3. Site-specific microbiome standards can be used to compare shotgun metagenomics and 16S rRNA analyses. Genomic DNA from the Vaginal Microbiome Whole Cell Mix (ATCC® MSA-2007™) was extracted and sequenced via shotgun whole-genome sequencing and 16S rRNA methods. The column on the left represents the theoretical (expected) abundance of each of the six organisms. The subsequent columns represent the actual number of reads observed per organism for each sequencing method. Here, two regions within the 16S rRNA gene were analyzed separately via two different Illumina® sequencing platforms (MiniSeq™ and MiSeq™).

Short-read versus Long-read Sequencing Platforms




Figure 4. Site-specific microbiome standards can be used to compare sequencing platforms. The Gut Microbiome Genomic DNA Mix (ATCC® MSA-1006™) was used to compare short-read and long-read sequencing platforms. The analysis of short-read sequencing data obtained on the Illumina MiSeq platform was consistent with the expected values. In contrast, the relative abundance varied among the long-read sequencing data obtained via the Oxford Nanopore GridION platform.

Data Analysis

Sequencing data was analyzed via the ATCC Microbiome Reference Standards module available through the One Codex data analysis platform. This module includes the following features:

- Accurate identification of reference organisms via local alignment analyses
- Conversion of read counts to relative abundance values by using genome size or 16S rRNA copy number
- Detection of non-reference organisms (false-positives) via the One Codex Database (whole-genome sequencing) or the Targeted Loci Database (16S rRNA)

Results are provided on a scorecard as true-positive, false-positive, and relative abundance scores.

Summary

This proof-of-concept study demonstrates the utility of site-specific microbiome standards as controls for evaluating run-to-run variability and optimizing assay performance at each stage of the microbiome analysis workflow.

- Whole cell standards can help identify biases introduced during DNA extraction and can also be used as full-process controls.
- Genomic DNA standards can be used for comparing various library preparation and sequencing platforms.
- The One Codex data analysis platform can be used to evaluate the number of true-positive, relative abundance, and false-positive scores for 16S rRNA community profiling and shotgun metagenomic sequencing methods.

Acknowledgements: We would like to thank Nick Greenfield, One Codex, for data analysis, Nezar Al-hebshi, Temple University, for oral microbiome design and Amanda Young and Jonathan Hetzel, Illumina, for their scientific input and sequencing expertise.

ATCC 10801 University Boulevard, Manassas, Virginia 20110-2209 **Phone:** 800.638.6597 **Email:** SalesRep@atcc.org **Web:** www.atcc.org

© 2019 American Type Culture Collection. The ATCC trademark and trade name, and any other trademarks listed in this publication are trademarks owned by the American Type Culture Collection unless indicated otherwise. Oxford Nanopore Technologies® (ONT) and GridION are trademarks of Oxford Nanopore Technologies. MiniSeq and MiSeq are registered trademarks of Illumina, Inc.