


# UNLOCK THE POTENTIAL OF AAV

AN INTRODUCTORY GUIDE TO  
ADENO-ASSOCIATED VIRUS (AAV) GENE DELIVERY



# TABLE OF CONTENTS

01

## Understanding AAV & Its Significance in Gene Delivery

- What is Adeno-associated Virus (AAV)?
- How Does AAV Compare to Other Viral Vectors?
- Benefits of AAV for Gene Delivery
- AAV Structure & Gene Delivery Mechanism

02

## Key Considerations in the Laboratory

- Producing AAV Particles in the Laboratory
- What is Multiplicity of Infection (MOI) & Why Does it Matter?
- AAV Biosafety

03

## Serotypes & Tissue Tropism

- Brief Overview
- Selecting Your Serotype

04

## AAV Packaging & Titration

- What Factors Effect Packaging Efficiency?
- Quantification of AAV Particles

05

## AAV Particle Storage

- What's the Best Way to Store AAV?

06

## Exploring OriGene's AAV Solutions:

- Products and Services Overview

Scan this QR code with your smartphone to view our AAV product portfolio




# Understanding AAV & Its Significance in Gene Delivery

## What is AAV?

Adeno Associated virus (AAV) is a small, non-enveloped virus belonging to the Parvoviridae family, which was first identified in 1965 as a contaminant of human and simian adenovirus preparations<sup>1</sup>. AAV has a single-stranded DNA genome of about 4.7kb in length, flanked by two inverted terminal repeats (ITRs) that serve as origins of replication and packaging signals.

This virus has gained significant attention in the research community due to its ability to deliver genetic material into targeted cells, with minimal immunogenicity and a broad range of infectivity. AAV remains a promising gene delivery vehicle for preclinical and clinical research and in recent years has made advances in a variety of therapeutic applications, such as treating Duchenne's Muscular Dystrophy (DMD).



*This virus has gained significant attention in the research community due to its ability to deliver genetic material into targeted cells, with minimal immunogenicity and a broad range of infectivity.*

## How Does AAV Compare to Other Viral Vectors?

Adeno-associated virus (AAV) stands out from other viral vectors like lentivirus and adenovirus in several ways. AAV exhibits a low immune response, enables long-term transgene expression, and carries a low risk of insertional mutagenesis. To delve deeper into the comparisons and contrasts among popular viral vectors used for gene delivery, please refer to the table below.

| Virus        | Expression           | Packaging Capacity | Genome         | Virus Size (nm) | Cells Infected            | Target Cell Genome Integration | Immune Response |
|--------------|----------------------|--------------------|----------------|-----------------|---------------------------|--------------------------------|-----------------|
| Lentivirus   | Stable               | <8 Kb              | RNA            | 80-130          | Dividing/<br>Non-dividing | Yes                            | Low             |
| AAV          | Transient or Stable* | ~4.5 Kb            | ssDNA (linear) | 18-26           | Dividing/<br>Non-dividing | No*                            | Very Low        |
| Adenovirus   | Transient            | >8 Kb              | dsDNA (linear) | 105             | Dividing/<br>Non-dividing | No                             | High            |
| y-Retrovirus | Stable               | <8 Kb              | RNA            | 80-130          | Dividing                  | Yes                            | Moderate        |

\*Recombinant AAV has a low frequency of target cell genome integration

Understanding the characteristics and differences among viral vectors used for gene delivery is crucial in designing effective gene therapy strategies. While each vector has its own advantages and considerations, AAV emerges as a particularly promising choice. By leveraging the unique properties of AAV and other viral vectors, researchers and clinicians can continue to advance the field of gene therapy, bringing us closer to transformative treatments for a wide range of genetic disorders.

Pre-made AAV Particles Available for 30,000+ Genes  
Human/Mouse/Rat

Browse Collection



## Benefits of AAV for Gene Delivery

Adeno-associated virus (AAV) is a promising viral gene delivery method with several key benefits that make it attractive for clinical and research applications. These benefits include:

- **Low immune response:** AAV naturally elicits a mild immune response in the body, which reduces the risk of host rejection or damage. This property of AAV can potentially prolong treatment time and limit the chance of unwanted side effects.
- **Long-term transgene expression:** AAV remains in the episome, therefore its expression is highly dependent on the cell's turnover rate, as this dilutes expression over time. The average time for expression is estimated to be 6-12 months<sup>2</sup>. This characteristic of AAV reduces the frequency of treatments required, providing a long-lasting therapeutic effect.
- **Broad range of infectivity:** Similar to other viral gene delivery methods, AAV can infect both dividing and non-dividing cells. This feature broadens its potential applications, as AAV can target a wide range of cell types in various diseases.
- **Low risk of insertional mutagenesis:** AAV does not integrate into the host genome, eliminating the potential risk of insertional mutagenesis. This significantly reduces the chances of unwanted complications associated with gene therapy treatments.

An illustrative case study highlighting the use of AAV in clinical settings is the development of Sarepta Therapeutics' gene therapy for Duchenne Muscular Dystrophy (DMD). On June 22nd, 2023, The FDA granted accelerated approval for ELEVIDYST<sup>™</sup>, the first gene therapy for Duchenne Muscular Dystrophy (DMD), which utilizes Adeno-associated Virus (AAV) to deliver mini-dystrophin gene to muscle tissues. This AAV-based therapy represents a significant advancement in gene therapy, showcasing the potential of AAV in treating genetic disorders.

## AAV Structure and Gene Delivery Mechanism

Recombinant AAV (rAAV) differs from wildtype-AAV in that it lacks viral DNA. rAAV has been described to act as a protein-based nanoparticle since it is optimized and engineered to traverse the cellular membrane to deliver its DNA to the nucleus<sup>3</sup>.

rAAV particles have the ability to form circular concatemers which persist as episomal DNA. Since this episomal DNA does not integrate into host genomes, it has the benefit of maintaining transient expression while avoiding insertional mutagenesis. This makes AAV stand out compared to other popular gene delivery methods such as lentivirus.

The process of AAV particle uptake through endocytosis is currently not fully understood. Among the different AAV serotypes, AAV2 has been extensively studied, providing the most available data on uptake mechanisms. Initially, research indicated that clathrin-mediated endocytosis is involved in the internalization of rAAV particles.

However, other studies have also shown evidence supporting alternative pathways that are independent of clathrin, as well as pathways involving GPI-enriched compartments. The data on endocytic mechanisms varies across different cell types, suggesting that these mechanisms may be dependent on the specific type of cell<sup>4</sup>.

Upon entry into the cell (Fig. 1), it is believed that AAV is transported through the Rab5+ early endosomal compartment, as well as the Rab7+ late endosome and Rab11+ recycling endosomes, although this can vary depending on the AAV serotype. After traversing the endomembrane system, AAV undergoes endosomal escape and remains in the cytosol.

The method by which AAV gains entry into the nucleus is still unknown. The variations observed in cell types and serotypes pose a significant challenge in comprehending the transduction mechanisms of AAV.

Once inside the nucleus, AAV undergoes uncoating, and its single-stranded DNA is converted into double-stranded DNA, which then undergoes transcription.




Fig 1. AAV Infection Mechanism

- References :
1. Rose JA, Hoggan MD, Shatkin AJ. Nucleic acid from an adeno-associated virus: chemical and physical studies. *Proc Natl Acad Sci U S A*. 1966 Jul;56(1):86-92. doi: 10.1073/pnas.56.1.86. PMID: 5229859; PMCID: PMC285679.
  2. Berns KI, Muzyczka N. AAV: An Overview of Unanswered Questions. *Hum Gene Ther*. 2017 Apr;28(4):308-313. doi: 10.1089/hum.2017.048. PMID: 28335618; PMCID: PMC5399733.
  3. Naso MF, Tomkowicz B, Perry WL 3rd, Strohl WR. Adeno-Associated Virus (AAV) as a Vector for Gene Therapy. *BioDrugs*. 2017 Aug;31(4):317-334. doi: 10.1007/s40259-017-0234-5. PMID: 28669112; PMCID: PMC548848.
  4. Berry GE, Asokan A. Cellular transduction mechanisms of adeno-associated viral vectors. *Curr Opin Virol*. 2016 Dec;21:54-60. doi: 10.1016/j.coviro.2016.08.001. Epub 2016 Aug 18. PMID: 27544821; PMCID: PMC5138113



# Key Considerations in the Laboratory

## Producing AAV Particles in the Laboratory

AAV viral particles are generated by co-transfecting typically a HEK293T cell line with a recombinant AAV vector carrying the gene of interest (GOI), a Rep/Cap plasmid (encoding AAV Rep/Cap genes), and an adenovirus-derived helper plasmid providing genes (E4, E2A, VA) needed for replication (Fig. 3).

**Day 1:** Plate HEK293T cells

**Day 2:** Co-transfect rAAV transfer plasmid, helper plasmids and Rep/Cap plasmid

**Day 4-7:** Harvest and Purify AAV particles

AAV particles can then be stored for later use at either 4°C for short term (<2 weeks), -80°C for long term (1 year) or immediately titered and used to transduce target cell lines.

## What is MOI & Why Does it Matter?

MOI stands for Multiplicity Of Infection and is the number of AAV particles added per cell during infection. Determining the appropriate MOI is important because it helps determine the optimal viral dose required for efficient transduction of your target cells, ensuring sufficient gene delivery and transgene expression levels.

The optimal MOI for each cell line varies significantly. AAV MOI ranges from 10,000 to 500,000 depending on serotype and cell type. The appropriate amount of viruses needed for the infection is crucial to the experimental result. Thus, it is strongly recommended to infect your target cells with an AAV control particle of your desired serotype, in your preliminary study to determine the optimal MOI.

## AAV Biosafety

Despite being derived from a non-pathogenic, wild-type virus, it is crucial to emphasize the significance of implementing robust biosafety measures when working with AAV. Understanding and adhering to stringent biosafety protocols is imperative to ensure the safe handling, production, and utilization of AAV for gene delivery in research and clinical settings.

OriGene highly recommends that you treat all AAV particles as BSL-2 organisms and strictly follow published BSL-2 guidelines with proper waste decontamination. AAV carrying harmful or toxic genes (e.g. activated oncogenes) require high levels of protection.

For more information on biosafety levels, please read the [NIH Biosafety Guidelines](#).





















































Fig 2. AAV Particle Production Infographic

# Serotypes and Tissue Tropism

Compared to other viral gene delivery methods, AAV's serotype-dependent tropism is a unique feature. These serotypes enable AAV to preferentially transduce specific organs or tissues. To date, there have been 12 natural serotypes and more than 100 variants of AAV isolated and studied as gene delivery vehicles<sup>1</sup>.

AAV2 is the most commonly used serotype. The natural serotypes include AAV1, AAV2, AAV3, AAV4, AAV5, AAV6, AAV7, AAV8, AAV9, AAV10 (AAVrh.10), AAV11, and AAV12.

The table<sup>2</sup> below demonstrates a few AAV serotypes and their various tropisms:

| Target Tissue  | Adeno-associated virus (AAV) Serotype  |  |  |  |  |  |  |  |
|--|--|--|--|--|--|--|--|--|
| <br><b>CNS</b>      |  AAV1   |  AAV2   |  AAV4   |  AAV5   |  AAV7     |  AAV9   |  AAV-rh10 |  AAV11      |
| <br><b>Liver</b>    |  AAV2   |  AAV3   |  AAV5   |  AAV7   |  AAV8     |  AAV9   |  AAV-rh10 |  AAV11      |
| <br><b>Lung</b>     |  AAV4   |  AAV5   |  AAV6   |  AAV9   |  AAV-rh10 |  AAV12  |  |  |
| <br><b>Muscle</b> |  AAV1 |  AAV2 |  AAV6 |  AAV7 |  AAV8   |  AAV9 |  AAV12  |  |
| <br><b>Heart</b>  |  AAV1 |  AAV4 |  AAV6 |  AAV7 |  AAV8   |  AAV9 |  |  |
| <br><b>Retina</b> |  AAV1 |  AAV2 |  AAV4 |  AAV5 |  AAV7   |  AAV8 |  AAV9   |  AAV-rh10 |

## Selecting your Serotype

Selecting the appropriate serotype is of the utmost importance since transduction efficiencies can be highly variable. Analyzing the literature is a good way to narrow down to a few targets, followed by testing with a panel of control particles of various serotypes that contain reporter constructs (on a small scale). This allows for visual analysis of transduction efficiencies. Testing kits are commercially available.

### References:

1.Li, C., Samulski, R.J. Engineering adeno-associated virus vectors for gene therapy. Nat Rev Genet 21, 255–272 (2020). <https://doi.org/10.1038/s41576-019-0205-4>

2.Issa, S.S.; Shaimardanova, A.A.; Solovyeva, V.V.; Rizvanov, A.A. Various AAV Serotypes and Their Applications in Gene Therapy: An Overview. Cells 2023, 12, 785. <https://doi.org/10.3390/cells12050785>

# AAV Packaging & Titration

## What Factors Effect Packaging Efficiency?

While there are many factors that contribute to packaging efficiency of AAV, these are the most common factors:

- DNA quality:** Potential recombination events can result in lower packaging efficiency and therefore must be tested. Plasmid recombination can be tested by carrying out a diagnostic enzyme digestion and Sanger sequencing.
- Optimal ratio of plasmids:** Studies report conflicting evidence for the optimal ratio for co-transfecting the required plasmids to produce AAV. OriGene recommends a 1:1:1 ratio<sup>1</sup>. i.e. Rep/Cap plasmid: Helper plasmid: Transfer plasmid
- Fragment length:** The packaging limit for AAV is around 4.7 kb, larger inserts may lead to lower packaged viral titer.
- Health of HEK293T cells:** HEK293T cells usually lose packaging efficiency after many passages, therefore, cells should not be used after culturing for 1-2 months. Secondly, cells are happier when seeded the day before transfection, which results in higher transfection efficiency.
- Transfection Reagent:** Since the Rep/Cap plasmid, Helper plasmid, and Transfer plasmid need to be co-transfected into HEK293T cells to produce viral particles, transfection efficiency is critical for high titer production. A transfection reagent that results in high transfection efficiency in HEK293T is essential for packaging.




Fig 3. In-vitro transduction efficiency of HEK293T cells with AAVDJ-CMV-GFP (Cat # CV900002S)

*Since the packaging limit for AAV is around 4.7 kb, larger inserts may lead to lower packaged viral titer.*

## Quantification of AAV Particles

Determining the titer of adeno-associated virus (AAV) is important in the production and characterization of AAV particles. The titer refers to the concentration or number of viral particles per unit volume.

Titering generally comes in two forms: physical or infectious.

Physical titer refers to the concentration of viral particles, also known as virions, in a sample. It quantifies the total number of viral particles, regardless of their infectivity. Physical titer is usually measured by techniques such as electron microscopy, spectrophotometry, quantitative polymerase chain reaction (qPCR), or Droplet Digital PCR (ddPCR).

Although qPCR is currently the most popular method used to quantify AAV, ddPCR can measure DNA copies without having to create a standard curve. The standard procedure for obtaining the physical titer involves lysing the cell and using PCR primers to bind to the ITR regions of the viral genome in the stock<sup>2</sup>.

Most commercial suppliers provide qPCR-based AAV titer kits.

Infectious titer, on the other hand, specifically measures the concentration of infectious viral units in a sample. AAV differs from other viral gene delivery systems, like lentivirus, in that the infectious titer is not reported due to its variability between serotypes and target cell lines.



### References :

1. Grieger JC, Soltys SM, Samulski RJ. Production of Recombinant Adeno-associated Virus Vectors Using Suspension HEK293 Cells and Continuous Harvest of Vector From the Culture Media for GMP FIX and FLT1 Clinical Vector. Mol Ther. 2016 Feb;24(2):287-297. doi: 10.1038/mt.2015.187. Epub 2015 Oct 6. PMID: 26437810; PMCID: PMC4817810.
2. Sanmiguell J, Gao G, Vandenbergh LH. Quantitative and Digital Droplet-Based AAV Genome Titration. Methods Mol Biol. 2019;1950:51-83. doi: 10.1007/978-1-4939-9139-6\_4. PMID: 30783968; PMCID: PMC7377043.



# AAV Particle Storage

## What's the Best Way to Store AAV?

In this section, we explore the key factors that influence AAV storage and provide you valuable tips to maintain the integrity and potency of your AAV preparations. Whether you are a seasoned researcher, a clinician working on gene therapy applications, or an industry professional involved in AAV production, these storage tips are essential for maximizing the shelf life and efficacy of your viral stocks.



01

### Minimize Free-Thaw Cycles

- Minimize freeze-thaw cycles. While AAV is relatively stable and can tolerate multiple cycles with minimal loss of activity, it is best to avoid unnecessary freeze-thawing.

02

### Avoid Drying

- Prevent drying of AAV samples, as it can result in protein denaturation.

03

### Prevent Introduction of Air

- Refrain from introducing air into the sample through vortexing, blowing bubbles, or similar operations, as it can lead to protein denaturation.

04

### Avoid Exposure to Denaturing Agents

- Avoid exposure to environmental extremes such as pH, chelating agents like EDTA, high temperatures, organic solvents, protein denaturants, and strong detergents as they can denature AAV.

05

### Optimal Long Term Storage

- For long-term storage, aliquot AAV samples and freeze them at -80°C if the virus will not be used within 1-2 weeks. This helps maintain the viability and activity of the AAV over an extended period.

06

### Maintain Appropriate Salt Concentrations

- Do not dilute AAV into low salt solutions. Some AAV serotypes, such as AAV-2, may aggregate in low salt environments, resulting in non-infectious aggregates if they are large. Maintain appropriate salt concentrations to prevent aggregation.

07

### Use Suitable Containers to Prevent Adhesion

- Minimize AAV's exposure to regular plastics, especially hydrophobic ones such as polystyrene, which can cause loss of activity. Store thawed AAV in siliconized or low protein binding tubes and use matching pipette tips. Add 0.01%-0.1% Pluronic F-68 to the formulation buffer to reduce sticking on regular plastics.



# Exploring OriGene's AAV Solutions: Products & Services

OriGene is committed to supporting and empowering researchers in their pursuit of scientific breakthroughs. Within this section, we showcase OriGene's comprehensive range of AAV solutions tailored to meet the diverse needs of the scientific community. From cutting-edge AAV vectors to personalized production services, OriGene offers an array of resources designed to enhance your research capabilities.


## Why Should You Use OriGene's Pre-made AAV-ORF™ Particles?

Our pre-made particles offer an exceptional platform for gene study with a host of advantageous features. With availability for over 30,000 targets, researchers gain unprecedented access to a diverse array of genes across various species, facilitating comprehensive investigations. Moreover, these particles are packaged with the AAV-2 serotype, a highly characterized and widely utilized serotype for gene delivery (options to customize are available). Another key attribute of these particles is their *Myc-DDK tag*, which is small enough to facilitate purification and detection without significantly limiting the gene insertion capacity.

### Key Benefits

- **Particles Arrive Transduction-ready:** Save up to 4 weeks of work
- **Ultra-Purified:** Mitigates risk of immune responses due to contaminants
- **High Titters:**  $>10^{13}$  GC/mL titer guaranteed for efficient transduction
- **Custom options available:** Flexible titers, serotypes, tags, and more, allowing you cater the AAV particles to your unique experiment

Scan the QR code to browse our collection pre-made AAV particles




## Your Research, Your Rules: Custom AAV Packaging Service


Producing high quality, high titer AAV can be difficult and time consuming. Our highly experienced AAV team is ready to tackle complex cloning & packaging projects to save you time and effort (Fig. 4).

### The OriGene Advantage

- Wide Serotype Selection: AAV1, AAV2, AAV3, AAV5, AAV6, AAV7, AAV8, AAV9, AAV-DJ, AAV-PHP.eb, AAVrh.10
- High Titters: Guaranteed titers  $>10^{13}$  GC/mL
- Flexibility Tagging Options: Myc-DDK, GFP, & RFP
- Fast Turnaround: 2-3 weeks



AAV 2



AAV DJ

Fig 4. Internal QC data of custom AAV2 & AAV-DJ particles using a GFP reporter.

## AAV Vectors

| Name                                      | Application     | Volume | Catalog #                       |
|---|-----------------|--------|---------------------------------|
| pAAV-AC-Myc-DDK Gene Expression Vector    | Over Expression | 10ug   | <a href="#"><u>PS100089</u></a> |
| pAAV-AC-GFP Gene Expression Vector        | Over Expression | 10ug   | <a href="#"><u>PS100090</u></a> |
| pAAV-AV-RFP Gene Expression Vector        | Over Expression | 10ug   | <a href="#"><u>PS100091</u></a> |
| pAAV-EF1a-tGFP-WPRE Gene Expression Vecor | Over Expression | 10ug   | <a href="#"><u>PS100127</u></a> |
| pGFP-A-shAAV Gene Expression Vector       | Knockdown       | 5ug    | <a href="#"><u>TR30034</u></a>  |

## AAV Control Particles

| Serotype   | Reporter | Volume     | Catalog #                        |
|------------|----------|------------|----------------------------------|
| AAV-1      | GFP      | 50ul/100ul | <a href="#"><u>CV900009S</u></a> |
| AAV-2      | GFP      | 50ul/100ul | <a href="#"><u>CV900001S</u></a> |
| AAV-3      | GFP      | 50ul/100ul | <a href="#"><u>CV900005S</u></a> |
| AAV-5      | GFP      | 50ul/100ul | <a href="#"><u>CV900011S</u></a> |
| AAV-6      | GFP      | 50ul/100ul | <a href="#"><u>CV900006S</u></a> |
| AAV-7      | GFP      | 50ul/100ul | <a href="#"><u>CV900008S</u></a> |
| AAV-8      | GFP      | 50ul/100ul | <a href="#"><u>CV900004S</u></a> |
| AAV-9      | GFP      | 50ul/100ul | <a href="#"><u>CV900007S</u></a> |
| AAV-rh.10  | GFP      | 50ul/100ul | <a href="#"><u>CV900010S</u></a> |
| AAV-DJ     | GFP      | 50ul/100ul | <a href="#"><u>CV900002S</u></a> |
| AAV-DJ     | P53-GFP  | 50ul/100ul | <a href="#"><u>CV900003S</u></a> |
| AAV-PHP.eb | GFP      | 50ul/100ul | <a href="#"><u>CV900012S</u></a> |



Find out more at [\*\*www.origene.com/products/cdna-clones/aav-products\*\*](http://www.origene.com/products/cdna-clones/aav-products)

© 2023 OriGene Technologies, Inc. All rights reserved. ELEVIDYS™ is the trademark of Sarepta Therapeutics. AAV-ORF and all other trademarks are the property of OriGene and its subsidiaries unless otherwise specified.

OriGene Technologies Inc.  
9620 Medical Center Drive, Suite #200  
Rockville, MD 20850  
USA

Phone: 1-301-340-3188  
Fax: 301-340-9254  
E-mail: [productorders@origene.com](mailto:productorders@origene.com)  
Website: [www.origene.com](http://www.origene.com)

